

Stainless Jugs



- Straight sided jugs with side handles
- 304 or 316L stainless steel construction
- Most sizes are graduated



Part No.	Stainless Steel	Nominal Volume	Graduated	Top Diameter	Height
A204-500	304	500ml	Yes	75mm	135mm
A204-1000	304	1000ml	Yes	100mm	135mm
A204-2000	304	2000ml	Yes	135mm	170mm
A204-3000	304	3000ml	No	120mm	230mm
A204-5000	304	5000ml	No	140mm	265mm
A206-500	316L	500ml	Yes	75mm	135mm
A206-1000	316L	1000ml	Yes	100mm	135mm
A206-2000	316L	2000ml	Yes	135mm	170mm

Jugs 316L ss



- Crevice free interior
- Polished finish
- 316L stainless steel



Part No.	Stainless Steel	Nominal Volume	Graduated	Top Diameter	Height
A186-500	316L	500ml	No	105mm	100mm
A186-1000	316L	1000ml	No	120mm	140mm
A186-2000	316L	2000ml	Yes	160mm	175mm

Conical Jugs




- Made from 304 stainless steel
- Crevice free design
- Polished finish
- Graduated



Part No.	Stainless Steel	Nominal Volume	Graduated	Top Diameter	Height
A154-500	304	500ml	Yes	115mm	115mm
A154-1000	304	1000ml	Yes	150mm	140mm
A154-2000	304	2000ml	Yes	171mm	200mm

316 Stainless Jugs (316Lss)

Volume 500ml to 5000ml



New

The 316 Stainless Steel Jugs are very robust and have a heavy duty design. These GMP correct jugs have straight sides, a pouring lip and a fully welded handle.

Each jug is highly polished and has a crevice-free interior. Graduation marks are etched onto the interior wall.

Part No.	Nominal Volume	Graduated	Diameter	Height	Material
A406-500	500ml	Yes	87mm	122mm	316L ss
A406-1000	1000ml	Yes	105mm	161mm	316L ss
A406-2000	2000ml	Yes	150mm	171mm	316L ss
A406-3000	3000ml	Yes	147mm	223mm	316L ss
A406-5000	5000ml	Yes	175mm	245mm	316L ss