

Metal Detectable Volumetric Spoon

Accurate measuring spoons



Metal Detectable Volumetric Spoons allow small volumes to be accurately measured. The Volumetric Spoons are designed to be used straight from the bag. Suitable for use with powders and liquids.

- Moulded & packed in a Cleanroom for superior cleanliness
- Individually Bagged
- Made from FDA, EC 1935/2004 & EU food grade 10/2011 conforming PP
- BSE/TSE Free
- Full Batch Traceability

NB. Volumetric spoons are not certified measuring devices. User validation is required. If in doubt please contact us



Part No.	Flat Fill Volume	Bowl Diameter	Overall Length	Material	Quantity per Pack
M1030B-25	2.5ml	23mm	158mm	Blue PP	4

Metal Detectable Long Handled Spoon

Long Reach Spoon



Our Metal Detectable Long Handle Spoons are ideal when you need to reach into deep containers. The long handle and the bright colour help them to be quickly and easily seen if they are left behind in the powder.

The ergonomic handle is comfortable to use with the spoon slightly angled.

- Moulded & packed in a Cleanroom for superior cleanliness
- Individually Bagged
- Made from FDA, EC 1935/2004 & EU food grade 10/2011 conforming PP
- BSE/TSE Free
- Full Batch Traceability



Part No.	Flat Fill Volume	Handle Length	Material	Quantity per Pack
M1036B-20	20ml	285mm	Blue PP	4

Metal Detectable Spoon

Highly Visible Sampling Spoon



The Metal Detectable Spoon from Sampling Systems is ideal for sampling small amounts of powder and for dispensing small quantities.

The bright blue colour ensures that the spoons are highly visible and easy to spot if left in most powders or granules (in line with HACCP guidelines).

Each spoon is individually bagged and can be used straight from the bag.

- Moulded & packed in a Cleanroom for superior cleanliness
- Individually Bagged
- Made from FDA, EC 1935/2004 & EU 10/2011 food grade conforming Polypropylene (PP)
- BSE/TSE Free
- Full Batch Traceability
- Available from stock



Part No.	Nominal Spoon Volume	Colour	Overall Length	Material	Quantity per Pack
M1006B-10	10ml	Blue	185mm	Blue PP	4

Metal Detectable

Metal Detectable Polypropylene (PP) - BLUE

Part No.	Nominal Spoon Volume	Sterile	Pack Quantity
M1006B-10	10ml	No	4



Metal Detectable Straight Spoon - PP (Polypropylene)

Part No.	Nominal Spoon Volume	Sterile	Pack Quantity
M1036B-20	20ml	No	4

