

Spoons reusable



Stainless Spoons 304ss

- 304 stainless steel construction
- Crevice free
- Bowl & flat ends
- Satin finish



Part No.	Nominal Bowl Volume	Nominal Overall Length	Material
A504-1	2ml	155mm	304 ss
A504-2	10ml	210mm	304 ss
A504-3	2ml	300mm	304 ss
A504-4	10ml	305mm	304 ss
A504-5	10ml	500mm	304 ss

Stainless Spoons 316L ss

- 316L stainless steel construction
- Crevice free
- Bowl & flat ends
- Satin finish



Part No.	Nominal Bowl Volume	Nominal Overall Length	Material
A506-300	10ml	300mm	316L ss
A506-600	10ml	600mm	316L ss
A506-1000	10ml	1000mm	316L ss

Volumetric Spoons



- Crevice free
- Satin finish
- Supplied as a set
- Each spoon is engraved with volume



Part No.	Stainless Steel	Volume	Bowl Diameter	Overall Length
A335-1	304	1.25ml	23mm	107mm
		5ml	33mm	112mm
		15ml	46mm	118mm
		100ml	79mm	127mm

Note: Volumetric Spoons are not certified measuring devices

Dosing and Volumetric Spoon

Volume 1.25ml to 100ml

New



Dosing and Volumetric Spoons are ideal for measuring powders and liquids. The range is available in a variety of sizes 1,25ml, 5ml, 15ml, 20ml and 100ml. The 5ml spoon is the same size as 1 teaspoon.

The spoons are made from high quality 304 stainless steel which means they have good chemical resistance.

Part No.	Volume	Bowl Diameter	Overall Length	Material
A335-012	1.25ml	23mm	106mm	304 ss
A335-050	5ml	33mm	112mm	304 ss
A335-150	15ml	47mm	119mm	304 ss
A335-200	20ml	51mm	141mm	304 ss
A335-1000	100ml	79mm	127mm	304 ss

Micro Spoon

New



- Ergonomic spoon
- Crevice free construction
- Made from 304 stainless steel

Part No.	Length	Bowl Width	Material
A634-150	150mm	5mm	304 ss
A634-200	210mm	5mm	304 ss

Open Spoon

New



- Single ended stainless steel spoon
- Pouring lip on end of bowl
- Crevice free construction
- Made from 304 stainless steel

Part No.	Length	Max. Width	Material
A620-200	200mm	25mm	304 ss

Chemical Spoon

Volume 5ml to 35ml

New



- Single bowl
- Flat handle
- Crevice free construction
- Made from 304 stainless steel

Part No.	Length	Nominal Volume	Material
A625-150	150mm	5ml	304 ss
A625-200	200mm	10ml	304 ss
A625-280	280mm	25ml	304 ss
A625-320	320mm	35ml	304 ss

Laboratory Spoon

Volume 2ml to 12ml

New



- Crevice free construction
- Made from 304 stainless steel

Part No.	Length	Nominal Volume	Material
A611-150	150mm	2ml	304 ss
A611-210	212mm	3ml	304 ss
A611-250	253mm	6ml	304 ss
A611-300	303mm	12ml	304 ss

Double Ended Laboratory Spoon

Volume 1ml to 12ml

New



- Double ended stainless steel spoon (different volume at each end)
- Flat handle
- Crevice free construction
- Made from 304 stainless steel

Part No.	Length	Nominal Volume	Material
A612-150	150mm	1ml & 2ml	304 ss
A612-210	212mm	1ml & 3ml	304 ss
A612-250	253mm	4ml & 7ml	304 ss
A612-300	303mm	7ml & 12ml	304 ss

Double Ended Spoon

Volume 0.5ml to 8ml

New



- Double ended stainless steel spoon (different volume at each end)
- Round handle
- Crevice free construction
- Made from stainless steel

Part No.	Length	Nominal Volumes	Material
A104-150	150mm	0.5ml & 2ml	Stainless steel
A104-180	180mm	0.8ml & 2.8ml	Stainless steel
A104-210	210mm	1.6ml & 3ml	Stainless steel
A104-250	250mm	1.8ml & 5ml	Stainless steel
A104-300	300mm	2ml & 8ml	Stainless steel