

Metal Detectable Liquid Funnel

Safely transfer liquids and suspensions

Our Metal Detectable Liquid Funnels are ideal when decanting liquids and suspensions. The Funnels have a small side handle to keep hands away from the product contact area when the funnels are in use.

- Moulded & packed in a Cleanroom for superior cleanliness
- Individually Bagged
- Made from FDA, EC 1935/2004 & EU 10/2011 conforming PP
- BSE/TSE Free
- Full Batch Traceability



Part No.	Top Diameter	Spout Diameter	Overall Height	Material	Quantity per Pack
M1023B-100	100mm	10mm	105mm	Blue PP	4