

Metal Detectable Powder Funnel

Wide spout for efficient powder flow



The Metal Detectable Powder Funnel is ideal for transferring powders and granules. The wide spout helps to prevent powder blockages. Metal Detectable Power Funnels are made from chemically resistant food grade PP.

- Individually bagged
- Made from FDA, EC 1935/2004 & EU 10/2011 conforming PP
- Moulded & packed in a Cleanroom for superior cleanliness
- Available from stock for fast delivery worldwide



| Part No. | Top Diameter | Spout Diameter | Material | Quantity per Pack |
|------------|--------------|----------------|----------|-------------------|
| M1051B-200 | 200mm | 43mm | Blue PP | 4 |