



Metal Detectable Ladle

Efficient surface sampling of liquids

The Metal Detectable Ladle is designed to quickly and efficiently take surface samples from liquids. They are moulded and packed in a Cleanroom and so can be used straight from the bag.

The one piece moulded design has been specially made to fit through a standard 2" lid bung. 2" bungs are standard features on many barrels and kegs.



- Moulded & packed in a Cleanroom for superior cleanliness
- Individually Bagged
- Made from FDA, EC 1935/2004 & EU 10/2011 conforming PP
- BSE/TSE free
- Full Batch Traceability
- Available from stock



Metal Detectable

Part No.	Nominal Volume	Diameter	Material of Construction	Quantity per Pack
M1033B-50	50ml	48mm	Blue PP	4