

Pallet Knives and Scrapers



Stirring Paddles


- 316L stainless steel
- Crevice free construction
- Satin finish
- Now available with holes for improved stirring



Part No.	Overall Length	Blade Width	Blade Length	Blade Type
A344-1000	1000mm	125mm	240mm	Plain
A344-1500	1500mm	150mm	300mm	Plain
A916-1000	1000mm	125mm	240mm	With Holes
A916-1500	1500mm	150mm	300mm	With Holes

GMP Pallet Knives

- Crevice free, all stainless steel construction
- Satin finish



Part No.	Blade Length	Blade Width	Overall Length
A330-100	100mm	25mm	185mm
A330-130	130mm	20mm	225mm
A330-150	150mm	25mm	255mm
A330-200	200mm	32mm	325mm
A330-250	250mm	38mm	380mm
A330-300	300mm	45mm	440mm
A330-350	350mm	50mm	495mm

Pallet Knife with Plastic Handle

- Stainless steel blade with plastic handle
- Satin finish



Part No.	Blade Length	Blade Width	Overall Length
A200-100	100mm	20mm	190mm
A200-200	200mm	30mm	310mm
A200-300	300mm	35mm	415mm

Chattaway Spatula

New



- Ergonomic spatula
- Crevice free construction
- Made from 304 stainless steel

Part No.	Length	Max. Width	Material
A627-100	100mm	9mm	304 ss
A627-150	150mm	9mm	304 ss
A627-210	210mm	9mm	304 ss
A627-250	250mm	10mm	304 ss
A627-300	303mm	10mm	304 ss

Micro Chattawat Spatula

New



- Ergonomic spatula
- Crevice free construction
- Made from 304 stainless steel

Part No.	Length	Max. Width	Material
A628-100	103mm	4mm	304 ss
A628-150	150mm	4mm	304 ss
A628-200	200mm	4mm	304 ss

Narrow Spoon Spatula

New



- Ergonomic spatula
- Crevice free construction
- Made from 304 stainless steel

Part No.	Length	Nominal Volume	Material
A618-235	235mm	2ml	304 ss

Curved Spatula

New



- Double ended spatula
- Crevice free construction
- Made from 304 stainless steel

Part No.	Length	Blade Width	Material
A630-120	120mm	15mm	304 ss
A630-210	210mm	25mm	304 ss

Pallet Knife with Wooden Handle



- Stainless steel blade with wooden handle
- Satin finish

Part No.	Blade Length	Blade Width	Overall Length
A201-100	100mm	20mm	190mm
A201-150	150mm	23mm	252mm
A201-200	200mm	30mm	310mm
A201-300	300mm	35mm	415mm

Stainless Steel Scrapers

- Crevice free, stainless steel construction
- Satin finish



Part No.	Stainless Steel	Blade Width	Overall Length
A350-75	420S45	75mm	232mm
A350-100	420S45	100mm	277mm
A656-100	316L	100mm	235mm

Flat Spatula

New



- Ergonomic spatula
- Crevice free construction
- Made from 304 stainless steel

Part No.	Length	Max. Width	Material
A631-120	120mm	21mm	304 ss
A631-210	210mm	25mm	304 ss

Flat Laboratory Spatula

New



- Double ended stainless steel flat spatula
- Crevice free construction
- Made from 304 stainless steel

Part No.	Length	Max. Width	Material
A613-150	149mm	20mm	304 ss
A613-210	212mm	25mm	304 ss
A611-250	255mm	32mm	304 ss
A611-300	306mm	40mm	304 ss