



## Metal Detectable Micro-Spatula

Ideal for handling small quantities of powders & granules



The Metal Detectable Micro-Spatula is ideal for sampling small volumes. Its narrow sampling blade enables samples to be taken from very small containers such as vials.

Alternatively it can be used for dispensing small amounts of powders and granules.

- Moulded & packed in a Cleanroom for superior cleanliness
- Individually Bagged
- Made from FDA, EC 1935/2004 & EU 10/2011 conforming PP
- BSE/TSE Free
- Full Batch Traceability



Part No.	Nominal Sample Volume	Material	Blade Width x Length	Overall Length	Quantity per Pack
M1013B-170	2ml	Blue PP	13 x 40mm	170mm	4

Metal Detectable

## Metal Detectable V-Spatula

Powder sampling made easy

Quick surface and core sampling is easy with our V-Spatula. Available in 2 lengths to offer shallow or deep sampling. They can also be used to take samples from sacks.

- Moulded & packed in a Cleanroom for superior cleanliness
- Ideal for general sampling and dispensing
- Suitable for free flowing and cohesive powders
- Ergonomic handle
- Free standing design



Part No.	Blade Length	Blade Width	Overall Length	Material	Quantity per Pack
M1010B-150	150mm	20mm	248mm	Blue PP	4
M1010B-250	250mm	20mm	345mm	Blue PP	4



## Metal Detectable Pallet Knife & Stirrer

Ideal for Production and Laboratory Use

The Metal Detectable Pallet Knife is very effective when working with creams and ointments, it is also a very efficient stirrer. The ergonomic handle makes them comfortable to use.

- Moulded & packed in a Cleanroom for superior cleanliness
- Individually Bagged
- Made from FDA, EC 1935/2004 & food grade EU 10/2011 conforming PP
- BSE/TSE Free
- Full Batch Traceability



Part No.	Blade Length	Blade Width	Overall Length	Material	Quantity per Pack
M1016B-200	200mm	45mm	310mm	Blue PP	4

## Metal Detectable Scraper

Cleanroom manufactured



The perfect tool for removing material from flat surfaces including tables, work tops and general equipment.

- Moulded & packed in a Cleanroom for superior cleanliness
- Individually Bagged
- Made from FDA, EC 1935/2004 & food grade EU 10/2011 conforming PP
- BSE/TSE Free
- Full Batch Traceability



Part No.	Blade Length	Blade Width	Overall Length	Material	Quantity per Pack
M1019B-110	110mm	110mm	225mm	Blue PP	4