

Heavy Duty Buckets (316L stainless Steel)



Stainless Buckets are made from heavy gauge 316L stainless steel and so are very durable.

- Made from 316L stainless steel
- Crevice free interior
- Swing handle
- Flat base
- Not graduated
- Lids sold separately



Part No.	Nominal Volume	Max Diameter	Height	Material
A586-10	10 litres	300mm	233mm	316L ss
A586-12	12 litres	310mm	262mm	316L ss
A586-15	15 litres	330mm	305mm	316L ss
A586-10L	Lid for 10L Bucket (A586-10) - 316L ss			
A586-12L	Lid for 12L Bucket (A586-12) - 316L ss			
A586-15L	Lid for 15L Bucket (A586-15) - 316L ss			

Buckets



- Crevice free interior
- Graduated
- Supplied with skirt band
- Robust design
- 304 stainless steel
- Lids sold separately (Lids not available for 6L bucket)



Bucket Part No.	Nominal Volume	Top Diameter	Height	Material
A264-6	6 litres	280mm	190mm	304 ss
A264-10	10 litres	300mm	295mm	304 ss
A264-12	12.5 litres	330mm	315mm	304 ss
A264-15	15 litres	350mm	360mm	304 ss
A264-10L	Lid for 10L Bucket (A264-10) - 304ss			
A264-12L	Lid for 12L Bucket (A264-12) - 304ss			
A264-15L	Lid for 15L Bucket (A264-15) - 304ss			

Buckets with Spouts (316L Stainless Steel)



- Bucket with pouring spout
- Made from 316L stainless steel
- Crevice free interior
- Not graduated
- Swing handle
- Supplied with skirt band
- Drop on lid. Covers the spout completely. Lids sold separately



Bucket Part No.	Nominal Volume	Overall Diameter	Nominal Height	Material
A686-15	15 litres	310mm	333mm	316L ss
A686-20	20 litres	330mm	387mm	316L ss
A686-15L	Lid for 15L Bucket with Spout (A686-15) - 316L ss			
A686-20L	Lid for 20L Bucket with Spout (A686-20) - 316L ss			

Medicine Cups (Small Volume)



Our medicine cups are ideal for putting small amounts of tablets, granules or powder into. The flat base design means they are ideal for weigh scales.

- 304 stainless steel
- Crevice free
- Polished



Part No.	Nominal Volume	Diameter	Height
A764-30	30ml	50mm	25mm

Bowls



- Crevice free design
- Polished finish
- Range of sizes
- 304 Stainless Steel



Part No.	Nominal Volume	Diameter	Height	Material
A534-1000	1000ml	170mm	85mm	304 ss
A534-2000	2000ml	220mm	110mm	304 ss
A534-4000	4000ml	270mm	130mm	304 ss

Buckets (316Lss)

Volume 10L to 15L

New



- Made from 316L stainless steel
- Crevice free interior
- Swing handle
- Flat base
- Non graduated

Part No.	Nominal Volume	Top Diameter	Height	Material
A586-10	10L	300mm	233mm	316L ss
A586-12	12L	310mm	262mm	316L ss
A586-15	15L	330mm	305mm	316L ss

Lids

Suitable for use with A586 buckets only

Part No.	To suit Bucket Volume	Material
A586-10L	10L	316L ss
A586-12L	12L	316L ss
A586-15L	15L	316L ss

Buckets (316Lss)

Volume 12L

New



- Crevice free interior
- Graduated
- Supplied with skirt band
- Robust design
- 316L stainless steel

Part No.	Nominal Volume	Top Diameter	Height	Material
A826-12	12L	310mm	285mm	316L ss

Bucket with Spout (316Lss)

Volume 15L to 20L

New



- Bucket with pouring spout
- Made from 316L stainless steel
- Crevice free interior
- Non graduated
- Swing handle
- Supplied with skirt band

Part No.	Nominal Volume	Overall Diameter	Height	Material
A686-15	15L	310mm	333mm	316L ss
A686-20	20L	330mm	387mm	316L ss