


Micro-Spatulas

- Stainless Steel
- Range of sizes
- Crevice free
- Bright polish finish



* A100-20 and A106-20



Part No.	Material	Blade Width	Blade Length	Overall Length
A100-1	304	9mm	50mm	210mm
A100-2	304	5mm	45mm	210mm
A100-3	304	2.5mm	35mm	150mm
*A100-20	304	12mm	40mm	175mm
*A106-20	316L	12mm	40mm	175mm

Trulla Spatula

New



- Special curved spatula
- Crevice free construction
- Made from 304 stainless steel

Part No.	Length	Max. Width of Ends	Material
A614-175	177mm	23mm	304 ss

Nuffield Spatula

New



- Ergonomic spatula
- Crevice free construction
- Made from 304 stainless steel

Part No.	Length	Max. Width of Ends	Material
A615-140	147mm	13mm	304 ss