

Instrument Trays



The Sampling Systems range of stainless steel trays are ideal for carrying tools etc.

- Range of sizes
- Available from stock
- Crevice free design
- High quality polished finish



Part No.	Dimensions (L x W x H)	Material
A674-215	215 x 160 x 18mm	304 Stainless
A674-295	295 x 197 x 18mm	304 Stainless
A674-395	395 x 275 x 11mm	304 Stainless
A674-490	490 x 320 x 15mm	304 Stainless

Stainless Trays - General Purpose



The general purpose range of trays are ideal for everything from mixing formulations to carrying waste.

- Range of sizes
- Available from stock
- Crevice free design
- High quality polished finish
- Durable



Part No.	Dimensions (L x W x H)	Material
A684-223	223 x 155 x 43mm	304 Stainless
A684-240	240 x 160 x 50mm	304 Stainless
A684-310	310 x 197 x 47mm	304 Stainless
A684-340	340 x 210 x 62mm	304 Stainless

Tablet Counting Triangle (316L ss)



New sizes

- Quick and easy way to count round tablets
- Suitable for all sizes of round tablets
- 316L stainless steel or PS shatter-proof plastic
- Easily used by right handed and left handed people
- Etched with a guide to show the number of tablets in the tray



Re-usable

Part No.	Side length	Height	Sterile	Material	Quantity per Pack
A356-160	160mm	14mm	No	316L ss	1
A356-275	275mm	12mm	No	316L ss	1

Disposable

Part No.	Side length	Height	Sterile	Material	Box Quantity
8066A-1	200mm	12mm	No	PS	30
8066S-1S	200mm	12mm	Yes	PS	30

Shallow Trays



- All purpose shallow tray
- Crevice free
- 304 stainless steel construction
- Satin finish



Part No.	Stainless Steel	Length	Width	Depth
A224-220	304	220mm	150mm	10mm
A224-290	304	290mm	210mm	10mm

Deep Trays



- Multi-purpose deep tray
- Crevice free
- 304 or 316L stainless steel construction
- Satin finish



Part No.	Stainless Steel	Length	Width	Depth
A234-220	304	220mm	150mm	45mm
A234-355	304	355mm	250mm	55mm
A234-395	304	395mm	280mm	70mm
A234-475	304	475mm	315mm	90mm
A234-530	304	530mm	325mm	100mm
A236-220	316L	220mm	150mm	45mm

Deep Trays (316Lss)

New



Stainless steel Deep Trays are a durable and hygienic alternative to plastic trays. They are made from high quality 316L stainless steel and are designed for heavy duty use.

Deep Trays (316L ss) are used in pharmaceutical, chemical and food companies around the world.

- 316L stainless steel
- Highly polished
- Rounded corners for easy cleaning

Part No.	Length	Width	Height	Material
A236-250	252mm	177mm	50mm	316L ss
A236-350	380mm	280mm	54mm	316L ss
A236-400	434mm	327mm	66mm	316L ss
A236-500	535mm	330mm	95mm	316L ss