

Metal Detectable Polypropylene (PP) - BLUE

Part No.	Nominal Volume	Sterile	Overall Height	Outer Diameter	Pack Quantity
M1047B-250	250ml	No	88mm	76mm	4
M1047B-500	500ml	No	106mm	94mm	4



Metal Detectable Bottle

Individually Bagged Bottles Made in a Cleanroom

Metal Detectable Bottles are designed to meet the requirements of the pharmaceutical, food and other hygienic industries.

The bottles are moulded and packed in a Class 100,000 Cleanroom. They have a wide opening and so are ideal when working with powders.

- Powder and splash tight
- Moulded & packed in a Cleanroom for superior cleanliness
- Made from FDA, EC 1935/2004 & food grade EU 10/2011 conforming PP
- BSE/TSE Free
- Full Batch Traceability
- Individually bagged



Part No.	Nominal Volume	Overall Height	Outer Diameter	Material	Quantity per Pack
M1047B-250	250ml	88mm	76mm	Blue PP	4
M1074B-500	500ml	106mm	94mm	Blue PP	4