

NEW

Scoops

from Sampling Systems



SteriWare PharmaScoop®

The Industry Standard - High Quality, Single Use Powder Scoops

The Sampling Systems range of disposable PharmaScoops® have been specially designed for dispensaries and manufacturing. They are the ultimate in quality and affordability.

SteriWare® PharmaScoops® are used by pharmaceutical, food and industrial companies all over the world.

- Manufactured and packed in a Cleanroom
- Individually Bagged in ultra strong packing with easy opening design
- Made from FDA, EC 1935/2004 & EU 10/2011 Conforming Materials
- BSE/TSE Free
- Full Batch Traceability
- Available sterile (SAL 10⁻⁶) and non-sterile from stock



The GMP PharmaScoop

It is Good Manufacturing Practice to keep your hands out of the product as much as possible. In accordance with this we have positioned the ergonomic handles high up on the scoop.

The Biggest Choice - Make sure you have the right scoop



Don't forget your scoop!
A brightly coloured scoop can be easily seen.

The Correct Material and Colour

We have the widest choice of materials and colours available:

PS - great for most applications. Strong and tough

HDPE and PP - extra levels of chemical resistance

Coloured scoops give increased visibility.

Wide Range of Sizes and Shapes

We know that not just one shape fits all applications so we have a wide range of styles and sizes so you can have the best scoop for your process.

Chemical Resistance Chart

See page 213.
Many commonly used chemicals are listed in our easy to use reference chart.

What Polymers are used in SteriWare products?

Abbreviation	Name	Approximate Working Temperature
HDPE	High Density Polyethylene	-50°C to +110°C
LDPE	Low Density Polyethylene	-50°C to +85°C
PP	Polypropylene	-0°C to +130°C
PS	Polystyrene (HIPS - High Impact Polystyrene)	-20°C to +80°C

Notes:

1. In this brochure PS and HIPS are used interchangeably. HIPS is used for most SteriWare products.
2. Working temperatures are approximate values only.

White PharmaScoop PS

Best Seller

Volume 30ml to 2500ml

Part No.	Nominal Capacity	Sterile	Box Quantity
8045A-30	30ml	No	100
8045A-30S	30ml	Yes	100
8045A-60	60ml	No	100
8045A-60S	60ml	Yes	100
8045A-125	125ml	No	100
8045A-125S	125ml	Yes	100
8045A-250	250ml	No	85
8045A-250S	250ml	Yes	85
8045A-500	500ml	No	50
8045A-500S	500ml	Yes	50
8045A-1000	1000ml	No	25
8045A-1000S	1000ml	Yes	25
8045A-2500	2500ml	No	15
8045A-2500S	2500ml	Yes	15



- ✓ Made from Polystyrene (PS)
- ✓ FDA and EC conforming
- ✓ BSE/TSE free
- ✓ Individually Bagged
- ✓ Made and packed in a Cleanroom

PharmaScoop HDPE

Extra Chemical Resistance

Volume 30ml to 2500ml

Part No.	Nominal Capacity	Sterile	Box Quantity
8045H-30	30ml	No	100
8045H-60	60ml	No	100
8045H-60S	60ml	Yes	100
8045H-125	125ml	No	100
8045H-125S	125ml	Yes	100
8045H-250	250ml	No	85
8045H-250S	250ml	Yes	85
8045H-500	500ml	No	50
8045H-500S	500ml	Yes	50
8045H-1000	1000ml	No	25
8045H-1000S	1000ml	Yes	25
8045H-2500S	2500ml	Yes	15



- ✓ Made from HDPE
- ✓ FDA and EC conforming
- ✓ BSE/TSE free
- ✓ Individually Bagged
- ✓ Made and packed in a Cleanroom

PharmaScoop PP

Great Chemical Resistance

Volume 30ml to 2500ml

Part No.	Nominal Capacity	Sterile	Box Quantity
8045P-30	30ml	No	100
8045P-30S	30ml	Yes	100
8045P-60	60ml	No	100
8045P-60S	60ml	Yes	100
8045P-125	125ml	No	100
8045P-125S	125ml	Yes	100
8045P-250	250ml	No	85
8045P-250S	250ml	Yes	85
8045P-500	500ml	No	50
8045P-500S	500ml	Yes	50
8045P-1000	1000ml	No	25
8045P-1000S	1000ml	Yes	25
8045P-2500	2500ml	No	15
8045P-2500S	2500ml	Yes	15



- ✓ Made from Polypropylene (PP)
- ✓ FDA and EC conforming
- ✓ BSE/TSE free
- ✓ Individually Bagged
- ✓ Made and packed in a Cleanroom

Red PharmaScoop PS

Highly Visible

Volume 30ml to 2500ml

Part No.	Nominal Capacity	Sterile	Box Quantity
8045R-30	30ml	No	100
8045R-30S	30ml	Yes	100
8045R-60	60ml	No	100
8045R-60S	60ml	Yes	100
8045R-125	125ml	No	100
8045R-125S	125ml	Yes	100
8045R-250	250ml	No	85
8045R-250S	250ml	Yes	85
8045R-500	500ml	No	50
8045R-500S	500ml	Yes	50
8045R-1000	1000ml	No	25
8045R-1000S	1000ml	Yes	25
8045R-2500	2500ml	No	15
8045R-2500S	2500ml	Yes	15



- ✓ Made from Polystyrene (PS)
- ✓ FDA and EC conforming
- ✓ BSE/TSE free
- ✓ Individually Bagged
- ✓ Made and packed in a Cleanroom

Blue PharmaScoop PP

Highly Visible with High Chemical Resistance

Volume 60ml to 2500ml

Part No.	Nominal Capacity	Sterile	Box Quantity
8045B-60	60ml	No	100
8045B-60S	60ml	Yes	100
8045B-125	125ml	No	100
8045B-125S	125ml	Yes	100
8045B-250	250ml	No	85
8045B-250S	250ml	Yes	85
8045B-500	500ml	No	50
8045B-500S	500ml	Yes	50
8045B-1000	1000ml	No	25
8045B-1000S	1000ml	Yes	25
8045B-2500	2500ml	No	15
8045B-2500S	2500ml	Yes	15



- ✓ Made from Polypropylene (PP)
- ✓ FDA and EC conforming
- ✓ BSE/TSE free
- ✓ Individually Bagged
- ✓ Made and packed in a Cleanroom

ULTRA PharmaScoop

USP Class VI - Double Bagged

- ✓ Made from HDPE - USP Class VI
- ✓ USP Class 661, FDA and EC conforming
- ✓ BSE/TSE free
- ✓ Double Bagged
- ✓ Made and packed in a Cleanroom

Part No.	Nominal Capacity	Sterile	Box Quantity
9045H-125S	125ml	Yes	50
9045H-250S	250ml	Yes	25
9045H-1000S	1000ml	Yes	10

Other sizes available on request



GMP Volumetric Scoop PS

Ergonomic Flat Scoop

Part No.	Nominal Capacity	Sterile	Box Quantity
8286A-25	25ml	No	100
8286A-25S	25ml	Yes	100
8286A-50	50ml	No	100
8286A-50S	50ml	Yes	100
8286A-100	100ml	No	100
8286A-100S	100ml	Yes	100
8286A-150	150ml	No	100
8286A-150S	150ml	Yes	100
8286A-250	250ml	No	100
8286A-250S	250ml	Yes	100
8286A-500	500ml	No	40
8286A-500S	500ml	Yes	40
8286A-1000	1000ml	No	30
8286A-1000S	1000ml	Yes	30

Volume 25ml to 1000ml



- ✓ Free-standing scoop with flat base
- ✓ Made from Polystyrene (PS)
- ✓ FDA and EC conforming
- ✓ BSE/TSE free
- ✓ Individually Bagged
- ✓ Made and packed in a Cleanroom

GMP Sample Scoop With Lid PS

Secures the Sample in the Scoop

Part No.	Nominal Capacity	Sterile	Box Quantity
8916A-25	25ml	No	50
8916A-25S	25ml	Yes	50
8916A-50	50ml	No	50
8916A-50S	50ml	Yes	50
8916A-100	100ml	No	50
8916A-100S	100ml	Yes	50

Volume 25ml to 100ml



- ✓ Click-on see through lid allows you to see sample
- ✓ Made from Polystyrene (PS)
- ✓ FDA and EC conforming
- ✓ BSE/TSE free
- ✓ Individually Bagged
- ✓ Made and packed in a Cleanroom

ULTRA Volumetric Scoop

USP Class VI - Double Bagged

- ✓ Made from HDPE - USP Class VI
- ✓ USP Class 661, FDA and EC conforming
- ✓ BSE/TSE free
- ✓ Double Bagged
- ✓ Made and packed in a Cleanroom

Part No.	Nominal Capacity	Sterile	Box Quantity
9286H-25S	25ml	Yes	50
9286H-50S	50ml	Yes	50
9286H-100S	100ml	Yes	50
9286H-150S	150ml	Yes	50
9286H-250S	250ml	Yes	25
9286H-500S	500ml	Yes	20
9286H-1000S	1000ml	Yes	20

USP Class VI
Double Bagged

Volume 25ml to 1000ml

New





Metal Detectable Scoop

Volume 30ml to 2500ml

New



Part No.	Nominal Volume	Colour	Material	Pack Quantity
M1000B-30	30ml	Blue	PP	4
M1000B-60	60ml	Blue	PP	4
M1000B-125	125ml	Blue	PP	4
M1000B-250	250ml	Blue	PP	4
M1000B-500	500ml	Blue	PP	4
M1000B-1000	1000ml	Blue	PP	4
M1000B-2500	2500ml	Blue	PP	4

Introducing SteriWare ULTRA

SteriWare ULTRA is the new premium range from Sampling Systems. All of the products in the SteriWare ULTRA range are double bagged and are made from HDPE that conforms to USP Class VI and USP 661.

All SteriWare ULTRA products are Cleanroom manufactured, sterilised by gamma, have a 5 year shelf life and come with full documentation.

- *Double Bagged*
- *USP Class VI*
- *Available from Stock*



SteriWare[®]
Ultra

Contact us for further information & prices

SteriWare® Long Handled Scoop

Keep hands away from the product

SteriWare Long Handled Scoops are the ideal tool for reaching into deep containers. The extra long 300mm handle will help operators keep their hands out of the product.

The long handle makes the scoop more visible if it is left in the powder. The red version helps this even further.

Both non-sterile and sterile scoops are available from stock. They are manufactured in a Cleanroom to ensure cleanliness and reduce particulate contamination.

- Manufactured and packed in a Cleanroom
- Individually Bagged
- Made from FDA, EC 1935/2004 & EU 10/2011 conforming HDPE
- BSE/TSE Free
- Available in uncoloured natural HDPE or bright red HDPE
- Full Batch Traceability
- Available sterile (SAL 10⁻⁶) and non-sterile

Volume 25ml to 100ml



New sizes

HDPE - Natural

Part No.	Nominal Scoop Volume	Sterile	Box Quantity
8341H-25	25ml	No	50
8341H-25S	25ml	Yes	50
8341H-50	50ml	No	50
8341H-50S	50ml	Yes	50
8341H-100	100ml	No	50
8341H-100S	100ml	Yes	50



HDPE - Red

Part No.	Nominal Scoop Volume	Sterile	Box Quantity
8341R-50S	50ml	Yes	50
8341R-100S	100ml	Yes	50




New

Metal Detectable Long Handled Scoop

Volume 100ml

New

Part No.	Nominal Volume	Overall Length	Overall Width	Material	Pack Quantity
M1040B-100	100ml	390mm	50mm	Blue PP	4





Open PharmaScoops are a popular alternative to our best selling PharmaScoops (page 141). Open PharmaScoops are designed to be easy to clean due to their crevice free design. They have a broad ergonomic handle which makes them comfortable to use even when wearing gloves.

Huge stocks for fast worldwide delivery.

- Heavy duty construction
- 316L stainless steel
- Crevice free
- Bright polish
- Rounded corners for easy cleaning
- GMP correct



Part No.	Nominal Capacity	Width	Body Length	Overall Length	Material
A336-250	250ml	100mm	160mm	280mm	316L stainless steel
A336-500	500ml	100mm	195mm	310mm	316L stainless steel
A336-750	750ml	125mm	205mm	320mm	316L stainless steel
A336-1000	1000ml	130mm	270mm	390mm	316L stainless steel
A336-1500	1500ml	150mm	235mm	350mm	316L stainless steel
A336-2000	2000ml	155mm	290mm	410mm	316L stainless steel
A336-4000	4000ml	205mm	310mm	430mm	316L stainless steel

PharmaScoops® (316L ss)

Volume 10ml to 5000ml



PharmaScoops are high quality stainless steel scoops. They are ideal for dispensing and for sampling.

- High quality scoops
- Heavy duty construction
- 316L stainless steel
- Crevice free
- Bright polish
- Rounded corners for easy cleaning
- GMP correct



PharmaScoops are a worldwide best seller in the pharmaceutical, food and chemical industries.

Part No.	Nominal Capacity	Body Diameter	Body Length	Overall Length	Material
A316-10	10ml	25mm	45mm	75mm	316L ss
A316-50	50ml	38mm	70mm	110mm	316L ss
A316-100	100ml	50mm	100mm	180mm	316L ss
A316-250	250ml	75mm	120mm	235mm	316L ss
A316-500	500ml	95mm	140mm	240mm	316L ss
A316-750	750ml	95mm	210mm	330mm	316L ss
A316-1000	1000ml	110mm	235mm	340mm	316L ss
A316-1250	1250ml	125mm	195mm	305mm	316L ss
A316-1700	1700ml	135mm	230mm	355mm	316L ss
A316-2500	2500ml	185mm	215mm	320mm	316L ss
A316-5000	5000ml	220mm	230mm	355mm	316L ss

Scorpion PharmaScoops®

Volume 1500ml to 2000ml



- Heavy duty construction
- 316L stainless steel
- Hygienic crevice free construction
- Satin finish
- Rounded corners for easy cleaning

Part No.	Nominal Capacity	Body Width	Body Length	Overall Height
A346-1500	1500ml	155mm	260mm	130mm
A346-2000	2000ml	165mm	300mm	150mm

Long Handle PharmaScoops®

Volume 10ml to 1000ml



- Extra Long handle
- Crevice free
- 316L stainless steel construction
- Rounded corners for easy cleaning
- GMP correct



Part No.	Nominal Capacity	Body Diameter	Body Length	Handle Length
A396-10	10ml	25mm	45mm	300mm
A396-50	50ml	38mm	70mm	300mm
A396-100	100ml	50mm	100mm	300mm
A396-500	500ml	95mm	140mm	600mm
A396-1000	1000ml	110mm	235mm	600mm

Weighing PharmaScoops®

Volume 500ml



- Free-standing scoop, ideal for use with weighing scales
- 316L stainless steel
- Rounded corners
- Crevice free construction
- Satin finish
- Robust handle

Part No.	Nominal Capacity	Maximum Width	Body Length	Overall Length
A376-500	500ml	110mm	205mm	320mm

Square PharmaScoops®

Volume 1500ml



- Heavy duty construction
- Made from 316L stainless steel
- Crevice free
- Satin finish
- Rounded corners for easy cleaning
- GMP correct

Part No.	Nominal Capacity	Width	Body Length	Overall Length
A366-1	1500ml	160mm	190mm	315mm

GMP Scoop

High quality, heavy duty powder scoops with improved ergonomics

Volume 750ml to 5000ml



GMP PharmaScoops are suitable for the most demanding of applications. The fully radiused body ensures they are easy to clean and minimises the risk of cross contamination.

For extra heavy duty work the GMP PharmaScoop is available with an optional Top Handle to improve the ergonomics.

- Heavy duty construction
- 316L stainless steel
- Hygienic crevice free construction
- Bright Polish
- Optional Top Handle
- Rounded corners for easy cleaning



Fully Welded

The scoop is fully welded with all the critical welds polished smooth




Top Handle (Optional)

Gives extra stability and leverage when scooping large volumes



Main Handle

Smooth and fully welded to the body



Robust Construction

Made from 1.2mm thick, 316L stainless steel

Radius

Smooth rounded back makes cleaning and operation easier

Polish

High quality polish all over

Part No.	Nominal Capacity	Top Handle	Body Diameter	Body Length	Overall Length	Material
A606-750	750ml	No	98mm	180mm	333mm	316L ss
A606-1500	1500ml	No	125mm	210mm	365mm	316L ss
A606-3000	3000ml	No	160mm	270mm	420mm	316L ss
A606-5000	5000ml	No	195mm	350mm	500mm	316L ss
A776-3000	3000ml	Yes	160mm	270mm	420mm	316L ss
A776-5000	5000ml	Yes	195mm	350mm	500mm	316L ss

SteriWare PharmaScoops®



SteriWare PharmaScoops® are high quality, single use, disposable plastic scoops. All SteriWare PharmaScoops® are manufactured and packed in a Cleanroom. They are widely used in the pharmaceutical and food industries all over the world.

SteriWare PharmaScoops® are ideal for use with potent drugs as well as highly coloured or perfumed products.

- Manufactured and packed in a Cleanroom to ensure cleanliness
- Individually wrapped
- Made from tough, virgin PS, HDPE or PP (FDA acceptable, EU 10/2011 and EC 1935/2004)
- Available non-sterile or pre-sterilised (validated procedure using gamma irradiation)
- Available from stock

Filling Scoop (316L ss)

Volume 750ml



The Filling Scoop can be used to fill bags and bottles directly through the handle of the scoop. The scoop has a special rounded funnel shape to help the material flow through the handle. This Filling Scoop is suitable for use with most free-flowing powders and granules.

- Empty scoop through the handle
- Heavy duty construction
- 316L stainless steel
- Hygienic crevice free construction
- Bright polish
- Rounded corners for easy cleaning

Part No.	Nominal Capacity	Width	Body Length	Overall Length
A436-750	750ml	101mm	183mm	313mm

Scorpion Scoop Eco

Volume 1000ml



- 304 stainless steel
- Satin finish
- Suitable for general, light duty use

Part No.	Nominal Capacity	Body Width	Body Length	Overall Height
A434-1000	1000ml	120mm	232mm	200mm

PTFE Scoops

Volume 100ml to 500ml



- Crevice free
- PTFE construction (FDA acceptable and USP Class VI)
- Rounded corners for easy cleaning
- GMP correct

Part No.	Nominal Capacity	Body Diameter	Body Length	Overall Length
A333-100	100ml	65mm	95mm	195mm
A333-500	500ml	100mm	140mm	265mm

Industrial Scoops

Volume 100ml to 1000ml



- Light duty construction
- Stainless steel
- Fully welded construction

Part No.	Nominal Capacity	Body Width	Body Length	Overall Length
A643-100	100ml	65mm	117mm	205mm
A643-200	200ml	85mm	140mm	225mm
A643-300	300ml	90mm	170mm	275mm
A643-500	500ml	105mm	185mm	290mm
A643-1000	1000ml	125mm	200mm	302mm

Economy Scoops

Volume 120ml to 2000ml




- Light duty construction
- 304 stainless steel

Important Note - Economy Scoops are not crevice free and are therefore not suitable for GMP applications

Part No.	Nominal Capacity	Diameter	Overall Length
A304-120	120ml	50mm	195mm
A304-450	450ml	75mm	255mm
A304-1000	1000ml	120mm	315mm
A304-1500	1500ml	120mm	345mm
A304-2000	2000ml	160mm	365mm

Ice Scoops

Volume 100ml



- Perforated to allow excess water to flow away
- 304 stainless steel

Part No.	Nominal Capacity	Body Diameter	Body Length	Overall Length
A314-100	100ml	50mm	100mm	180mm

Available for Download
from www.sampling.com

New



The Range Brochure



Stainless Ware Brochure



FreeGlide Brochure



Metal Detectable Brochure



SteriWare Brochure



Label Solutions Brochure

 **Via Torre, 2**
 **28040 Borgo Ticino (NO) - Italy**
 **+39 0321 908466**
 **customercare@exonder.it**
 **<https://www.exonder.it>**



Sampling Systems Ltd.
4 & 5 Forge Mills Park,
Coleshill,
Warwickshire,
B46 1JH
UK

www.sampling.com
Tel: +44 (0)1675 466992
Email: info@sampling.com

QMS[®]
ISO 9001
REGISTERED

VISA



'SteriWare', 'PharmaScoop' and 'PharmaBottle' are registered Trade Marks belonging to Sampling Systems Ltd.

All goods are sold subject to our normal Terms & Conditions. Contact us for a copy.

©2023all rights reserved