


S.C.A.T. Europe Overview

Range of Canisters



Product Group 1: Electrically Conducting Canisters

Image	Thread	Capacity	Material	Colour	UN-Certified	Special Feature	Dimensions mm W x H x D
A	S 50	5 litres	PE-HD-el	black	No	designed as space-saving canister	65 x 335 x 335
B	S 50	10 litres	PE-HD-el	black	Yes	UN-X approval; with viewing strip	190 x 315 x 230
C	S 60 / 61	10 litres	PE-HD-el	black	Yes	UN-Y approval	185 x 265 x 290
D	S 60 / 61	10 litres	PE-HD-el	black	No	with red viewing strip	185 x 265 x 290
E	S 60 / 61	10 litres	PE-HD-el	black	No	with green viewing strip	185 x 265 x 290
F	S 60 / 61	10 litres	PE-HD-el	black	No	with yellow viewing strip	185 x 265 x 290
G	S 60 / 61	10 litres	PE-HD-el	black	No	with blue viewing strip	185 x 265 x 290
H	S 60 / 61	10 litres	PE-HD-el	black	No	with float	185 x 265 x 290
I	S 60 / 61	10 litres	PE-HD-el	black	No	with holder for sensor 108 178	185 x 265 x 290
J	S 60 / 61	20 litres	PE-HD-el	black	Yes	UN-Y approval	185 x 500 x 290
K	S 60 / 61	20 litres	PE-HD-el	black	No	with float	185 x 500 x 290
L	S 60 / 61	30 litres	PE-HD-el	black	Yes	UN-Y approval	240 x 455 x 364
M	S 60 / 61	30 litres	PE-HD-el	black	No	with float	240 x 455 x 364
N	S 70 / 71	60 litres	PE-HD-el	black	Yes	UN-Y approval	330 x 625 x 396
O	S 70 / 71	60 litres	PE-HD-el	black	No	with float	330 x 690 x 395
P	S 90	10 litres	PE-HD-el	black	Yes	UN-X approval; with viewing strip	195 x 380 x 195

S.C.A.T. Europe Overview

Range of Canisters



THE System for
Solvent Safety



Product Group 2: Non-Electrically Conducting Canisters


Image	Thread	Capacity	Material	Colour	UN-Certified	Special Feature	Dimensions mm W x H x D
A	GL 45	2,5 litres	PE-HD	natural (white)	Yes	UN-Y approval	115 x 210 x 150
B	GL 45	5 litres	PE-HD	natural (white)	Yes	UN-Y approval	150 x 250 x 195
C	GL 45	10 litres	PE-HD	natural (white)	Yes	UN-Y approval	190 x 305 x 230
D	S 51	5 litres	PE-HD	natural (white)	Yes	UN-X approval	145 x 250 x 190
E	S 55	5 litres	PE-HD	natural (white)	Yes	UN-X approval	182 x 240 x 162
F	S 55	5 litres	PE-HD	natural (white)	Yes	UN-Y approval	160 x 230 x 185
G	S 55	5 litres	F-PE-HD	natural (white)	Yes	UN-X approval; with inner fluorination	145 x 250 x 190
H	S 55	10 litres	F-PE-HD	natural (white)	Yes	UN-X approval; with inner fluorination	190 x 310 x 230
I	S 55	10 litres	PE-HD	natural (white)	Yes	UN-Y approval	185 x 305 x 225
J	S 60 / 61	5 litres	PE-HD	natural (white)	Yes	UN-X approval	164 x 235 x 185
K	S 60 / 61	12 litres	PE-HD	natural (white)	Yes	UN-X approval	200 x 350 x 235
L	S 60 / 61	20 litres	PE-HD	natural (white)	Yes	UN-X approval	260 x 390 x 289
M	S 60 / 61	20 litres	PE-HD	natural (white)	No	with float	260 x 455 x 285
N	S 60 / 61	30 litres	PE-HD	natural (white)	Yes	UN-Y approval	290 x 400 x 380
O	S 60 / 61 & GL 45	30 litres	PE-HD	natural (white)	No	2 openings	290 x 400 x 380
P	S 70 / 71	60 litres	PE-HD	natural (white)	Yes	UN-Y approval	330 x 635 x 370
Q	S 90	10 litres	PE-HD	natural (white)	Yes	UN-Y approval	195 x 380 x 195
R	S 95	5 litres	PE-HD	natural (white)	Yes	UN-Y approval	170 x 310 x 170
S	S 95	20 litres	PE-HD	natural (white)	No		260 x 390 x 290

S.C.A.T. Europe Overview

Range of Canisters



THE System for
Solvent Safety



107 334



107 335



107 335



107 336

Product Group 3: Space-Saving Canisters

Image	Thread	Capacity	Material	Colour	UN-Certified	Special Feature	Dimensions mm W x H x D
A	S 50	5 litres	PP	natural (white)	No		65 x 335 x 330
B	S 50	5 litres	PP	natural (white)	No	with float	65 x 330 x 405

Product Group 4: Politainers

Image	Thread	Capacity	Material	Colour	UN-Certified	Special Feature	Dimensions mm W x H x D
C	GL 38	5 litres	PE-LD / EVA Copolymer	natural (white)	Yes, but only together with Box 107 334	Box with UN-Y approval	178 x 178 x 178
D	GL 38	10 litres	PE-LD / EVA Copolymer	natural (white)	Yes, but only together with Box 107 335	Box with UN-Y approval	228 x 228 x 228
E	S 60 / 61	10 litres	PE-LD / EVA Copolymer	natural (white)	Yes, but only together with Box 107 335	Box with UN-Y approval	228 x 228 x 228
F	S 60 / 61	20 litres	PE-LD / EVA Copolymer	natural (white)	Yes, but only together with Box 107 336	Box with UN-Y approval	285 x 285 x 285

S.C.A.T. Europe Overview

Range of Canisters



THE System for
Solvent Safety



Product Group 5: Round Canisters

Image	Thread	Capacity	Material	Colour	UN-Certified	Special Feature	Dimensions mm W x H x D
A	B 53	2 litres	PE-HD	natural (white)	No		119 x 260 x 119
B	S 65	5 litres	PE-HD	natural (white)	No	with handle and measuring scale	167 x 330 x 167
C	S 65	10 litres	PE-HD	natural (white)	No	with handle and measuring scale	205 x 430 x 205
D	S 65	25 litres	PE-HD	natural (white)	No	w. carrying handles and measuring scale	278 x 580 x 278
-	S 65	60 litres	PE-HD	natural (white)	No	w. carrying handles and measuring scale	350 x 825 x 350
E	B 83	4 litres	PE-HD	natural (white)	No		155 x 338 x 155
F	B 83	10 litres	PP	natural (white)	No	with carrying handles	250 x 390 x 250
G	B 83	50 litres	PP	natural (white)	No	with carrying handles	380 x 680 x 380

Product Group 6: Flat Canisters

Image	Thread	Capacity	Material	Colour	UN-Certified	Special Feature	Dimensions mm W x H x D
H	B 63	8 litres	PP	natural (white)	No	with drain valve	375 x 165 x 305
I	S 70 / 71	20 litres	PE-HD	natural (white)	Yes	UN-Y approval	283 x 292 x 376