

# SCAT Europe Overview

## Range of Laboratory Glass Bottles



THE System for  
Solvent Safety



### Laboratory Glass Bottles

Image	Description	Thread	Capacity	Glass type	Form	Protective coating	Part Number
<b>A</b>	Original GL45 Screw Cap	GL 45	-	-	-	-	199 011
<b>B</b>	Laboratory bottle DURAN, GL45, 250 ml, Type 1	GL 45	250 ml	Clear glass	round	yes	101 994
<b>C</b>	Laboratory bottle DURAN, GL45, 250 ml, Type 3	GL 45	250 ml	Clear glass	round	-	501 117
<b>D</b>	Laboratory bottle DURAN, GL45, 250 ml, Type 2	GL 45	250 ml	Clear glass	square	-	501 112
<b>E</b>	Laboratory bottle DURAN, GL45, 500 ml, Type 1	GL 45	500 ml	Clear glass	round	yes	101 995
<b>F</b>	Laboratory bottle DURAN, GL45, 500 ml, Type 3	GL 45	500 ml	Clear glass	round	-	501 116
<b>G</b>	Laboratory bottle DURAN, GL45, 500 ml, Type 4	GL 45	500 ml	Brown glass	round	-	501 120
<b>H</b>	Laboratory bottle DURAN, GL45, 500 ml, Type 2	GL 45	500 ml	Clear glass	square	-	501 115
<b>I</b>	Laboratory bottle DURAN, GL45, 1 L, Type 1	GL 45	1 litre	Clear glass	round	yes	101 996
<b>J</b>	Laboratory bottle DURAN, GL45, 1 L, Type 2	GL 45	1 litre	Clear glass	round	-	501 113
<b>K</b>	Laboratory bottle DURAN, GL45, 1 L, Type 3	GL 45	1 litre	Brown glass	round	-	501 119
<b>L</b>	Laboratory bottle DURAN, GL45, 1 L	GL 45	1 litre	Clear glass	square	-	501 110

#200web / 01092021

# SCAT Europe Overview

## Range of Laboratory Glass Bottles



THE System for  
Solvent Safety



### Laboratory Glass Bottles

Image	Description	Thread	Capacity	Glass type	Form	Protective coating	Part Number
<b>M</b>	Laboratory bottle DURAN, GL45, 2 L, Type 1	GL 45	2 litres	Clear glass	round	yes	101 997
<b>N</b>	Laboratory bottle DURAN, GL45, 2 L, Type 2	GL 45	2 litres	Clear glass	round	-	501 118
<b>O</b>	Laboratory bottle DURAN, GL45, 5 L, Type 1	GL 45	5 litres	Clear glass	round	yes	101 998
-	Laboratory bottle DURAN, GL45, 5 L, Type 2	GL 45	5 litres	Clear glass	round	-	501 125
<b>P</b>	Laboratory bottle DURAN, GL45, 10 L	GL 45	10 litres	Clear glass	round	-	501 126
<b>Q</b>	Laboratory bottle DURAN, GLS80, 500 ml, Type 1	GLS80	500 ml	Clear glass	round	-	501 150
<b>R</b>	Laboratory bottle DURAN, GLS80, 1 L, Type 1	GLS80	1 litre	Clear glass	round	-	501 151
<b>S</b>	Laboratory bottle DURAN, GLS80, 1 L, Type 2	GLS80	1 litre	Brown glass	round	-	501 157
<b>T</b>	Laboratory bottle DURAN, GLS80, 2 L, Type 1	GLS80	2 litres	Clear glass	round	-	501 152

#200web / 01092021