




THE System
for Solvent Safety

The next Level!

SafetyCaps

Safe Solvent Supply for
any HPLC / UHPLC System



2.0



SafetyCaps 2.0

The next Level!



THE System
for Solvent Safety

NEW: Luer Lock Adapter
Useful for many applications

NEW: Air valve
for analytical **and**
preparative applications




NEW: PFA Fittings
Improved design for
easy tubing connection
Excellent chemical resistance
Flammability classification
V-0 (UL-94)

NEW: Medical
Grade PTFE
Excellent resistance
Food safe
FDA conform
Flammability
classification V-0
(UL-94)

NEW: Insertion Label
Mark your bottles with individual
descriptions



NEW: PPS Screw Cap
Autoclavable / sterilisable up to 200 °C
Flammability classification V-0 (UL-94)



* Air valve and Fittings for tubing connection
are included in delivery.



| | Description | Unit | Thread size | Connections for Ø 3.2 mm (1/8") OD tubing | Blind plugs for tubing connections |
|--|--|------|--------------|---|--|
| | SafetyCap I (V 2.0) * | 1 | GL 45 | 1x | - |
| | SafetyCap II (V 2.0) * | 1 | GL 45 | 2x | 1x |
| | SafetyCap III (V 2.0) * | 1 | GL 45 | 3x | 2x |
| | SafetyCap IV (V 2.0) * | 1 | GL 45 | 4x | 2x |
| | SafetyCap VI (V 2.0) * | 1 | GL 45 | 6x | 2x |
| | Air valve for SafetyCaps | 1 | UNF 1/4" 28G | - | - |
| | Air valve for SafetyCaps | 10 | UNF 1/4" 28G | - | - |
| | Fitting for Ø 3.2 mm OD tubing, blue | 5 | UNF 1/4" 28G | - | - |
| | Fitting for Ø 3.2 mm OD tubing, black | 5 | UNF 1/4" 28G | - | - |
| | Fitting for Ø 3.2 mm OD tubing, yellow | 5 | UNF 1/4" 28G | - | - |
| | Blind plug for tubing connections | 10 | UNF 1/4" 28G | - | - |



Phone

+49 (0) 61 05 / 30 55 86 - 0

e-Mail

info@scat-europe.com

Internet

www.scat-europe.com

SafetyCaps 2.0

HPLC Starter Kits



- Price advantage compared to single order
- Suitable for all common HPLC systems (independent of manufacturer)
- Use blind plugs (included) to lock unused connections



| Fig. | Part No. | Description | Content | Unit | Thread size | Connections for Ø 3.2 mm OD tubing |
|----------|------------------------------|-------------|----------------------------------|------|--------------|------------------------------------|
| A | Starter-Kit-Supply-S1 | | SafetyCap I (V 2.0) | 3x | GL 45 | 1 x |
| | | | SafetyCap II (V 2.0) | 1x | GL 45 | 2x |
| | | | Blind plug for tubing connection | 1x | UNF 1/4" 28G | - |
| B | Starter-Kit-Supply-S2 | | SafetyCap II (V 2.0) | 4x | GL 45 | 2x |
| | | | Blind plug for tubing connection | 4x | UNF 1/4" 28G | - |



THE System
for Solvent Safety

S.C.A.T. Europe GmbH
Waldecker Straße 7
64546 Mörfelden-Walldorf
Germany

☎ Phone: +49 (0) 61 05 / 30 55 86 - 0

✉ e-Mail: info@scat-europe.com

🌐 Internet: www.scat-europe.com

📍 Via Torre, 2
📍 28040 Borgo Ticino (NO) - Italy
☎ +39 0321 036125
✉ customercare@vivaxsrl.it
🌐 <https://www.vivaxlab.com>

