

# SafetyWasteCap 2.0

## S 60/61 Thread

### Universal



THE System  
for Solvent Safety

#### ORIGINAL SAFETY:

Innovative & universally usable

- SafetyWasteCap with S 60/61 thread in the new design of generation 2.0.
- More connections thanks to the beveled core - for more applications!
- Maximum flexibility and equipment on a single cap - no adapter required!
- Connections which are not required can be closed with blind plugs.
- Large scope of delivery:  
9x Fittings: 3x 1.6 mm, 3x 2.3 mm, 3x 3.2 mm  
3x Blind plugs for capillary connection  
3x Tube connectors  
3x Blind plugs for tube connections



#### NEW: Ergonomic Design

- Enhanced stability and construction / Improved handling / Easy installation / Smart exchange

#### NEW: Medical Grade PTFE

- Excellent chemical resistance / Food safe / FDA compliant / Flammability classification V-0 (UL-94)

#### NEW: PPS Screw Cap

- Autoclavable / sterilisable up to 200 °C / Flammability classification V-0 (UL-94)

#### Proven: PTFE Body with integrated FKM Seal

- Provides hermetical connection to the container opening and adapts smoothly even to rugged glass surfaces. An additional PTFE layer protects the seal from direct contact to solvents.



☎ Phone: +49 (0) 61 05 / 30 55 86 - 0

✉ e-Mail: [info@scat-europe.com](mailto:info@scat-europe.com)

🌐 Internet: [www.scat-europe.com](http://www.scat-europe.com)

S.C.A.T. Europe GmbH  
Waldecker Straße 7  
64546 Mörfelden-Walldorf  
Germany



Description	Thread size	Connections for Ø 1.6 / 2.3 / 3.2 mm OD tubing	Connections for Ø 6.4 - 9 mm ID tubing	Service life (in months)	Unit
SafetyWasteCap S 60/61, V2.0, Universal	S 60/61	3x	3x	-	1
Exhaust filter size M, V2.0, with splash protection and change indicator <b>Economy package</b>	GL 14	-	-	2x 6	2