

HPLC Starter Kit 3 V2.0

SafetyCap 2.0


GL 45 Thread

Universal!



THE System
for Solvent Safety

Initial equipment for your HPLC system -
Safe extraction of solvents.



- Price advantage compared to ordering piece by piece.
- Suitable for all current HPLC systems.
- Unused connectors can be closed by the blind screws delivered with the set.
This makes the system universally applicable, even when less connectors are needed for a time.



	Description	Thread size	Connections for Ø 1.6 / 2.3 / 3.2 mm OD tubing	Service life (in months)	Colour	Sales Unit ----- incl.
	HPLC Starter Kit 3, V2.0					1
incl.	SafetyCap III, V2.0	GL 45	3		blue/white	4x
incl.	Air valve for SafetyCaps, V2.0	UNF 1/4" 28G		6 months each	red	4x
incl.	Fitting 1.6 mm	UNF 1/4" 28G			green	12x
incl.	Fitting 2.3 mm	UNF 1/4" 28G			violet	12x
incl.	Fitting 3.2 mm	UNF 1/4" 28G			blue	12x
incl.	Blind plug	UNF 1/4" 28G			translucent	8x
incl.	Insertion Label / Title Block				white	8x